

Simvastatin Rosemont Oral Suspensions

20mg/5mL and 40mg/5mL

Instructions for administration via nasogastric (NG) or percutaneous endoscopic gastrostomy (PEG) tubes:
Simvastatin Rosemont Oral Suspension is suitable for use with the following types of NG and PEG tubes:

Material	External Bore Size (Fr Unit)	Internal Diameter (mm)	Maximum Length (cm)
Silicone	4	0.8	125
	6	1.0	125
	10	2.0	125
PVC	4	0.8	125
	8	1.5	125
	12	2.5	125
Polyurethane	4	0.8	125
	8	1.5	125
	12	2.6	125
	18	4.0	125

Ensure that the enteral feeding tube is free from obstruction before administration.

1. Flush the enteral tube with sufficient water to remove any residual feed.

2. Administer the required dose of Simvastatin Rosemont Oral Suspension with a suitable measuring device.

3. Flush the enteral tube with 10mL of water.

This product has not been tested with latex NG or PEG tubes and therefore should not be used for administration with tubes made from latex.

Air flushing of NG or PEG tubes to deliver the required dose of Simvastatin is not recommended.

Any unused product or waste material should be disposed of in accordance with local requirements.



DTM063.1 July 2020

Rosemont Pharmaceuticals Ltd. Rosemont House, Yorkdale Industrial Park, Braithwaite Street, Leeds LS11 9XE **T** +44 (0)113 244 1400 **F** +44 (0)113 245 3567 **E** infodesk@rosemontpharma.com Sales/ Customer Service: **T** +44 (0) 113 244 1999 **F** +44 (0) 113 246 0738 **W** www.rosemontpharma.com

Adverse events should be reported. Reporting forms and information can be found at www.mhra.gov.uk/yellowcard Adverse events should also be reported to Rosemont Pharmaceuticals Ltd on 0113 244 1400.